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 **Instruction Manual**

DDBioSplint C

DDBioSplint C

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Material-related characteristics

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Polycarbonat (PC) | 100

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Biegefestigkeit	[MPa]	≥ 98
Biegemodul	[MPa]	> 2300
Wasseraufnahme	[µg/mm ³]	3,6
Löslichkeit	[µg/mm ³]	0,2
Restmonomergehalt	[%]	≤ 0,0015

Composition [wt.%]

Polycarbonate (PC) | 100

Physical specifications

Flexural strength	[MPa]	≥ 98
Flexural modulus	[MPa]	> 2300
Water absorption	[µg/mm ³]	3.6
Solubility	[µg/mm ³]	0.2
Residual monomer content	[%]	≤ 0.0015

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1. Indications

DD Bio Splint C are transparent dental milling blanks out of polycarbonate for the manufacture of splints, therapeutic splints, bite regulators and drilling templates for short-term application in the oral cavity up to 30 days.

2. Contraindications

- intolerance to the contained components
- bridges with two molar pontics
- too little connector cross section
- model casting with clasps, bars, narrow sublingual bars and transversal connectors

3. Safety information

Please pay attention to the information in the instruction manual in its current version. Avoid the inhalation of dust particles during processing. Wear protection gloves, safety glasses and a dust mask to prevent irritations of eyes, skin and respiratory system.

4. Handling and storage

Only use clean, defect-free material out of undamaged packaging and make sure that the blanks are always stored in their original container in a dry and cool place, protected from light. Avoid contact with fluids.

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Instruction Manual

5. Processing

Handling this medical device should only be allowed to trained staff.
In general all constructions should be designed in a way to enable the compensation and the transfer of the mastication forces to the abutment teeth.

6. Milling and further processing

Only use one-edged milling tools in the milling machine or wet grinding procedures!

Further processing:

Tungsten carbide burs like the one you use for soft relin materials are suitable for rough adaption. For the finalization common cross-toothed millers can be used. Polishing should be carried out with regard to a low contact pressure, otherwise the resulting heat might lead to smearing or melting.

7. Cleaning

- ultrasonic bath (max. 40 °C) with mild cleaning agent
- process of sterilisation in an autoclave up to 121 °C

Do not use the following methods for cleaning:

- steam blasting
- organic solvents and acids (e.g. acetone, hydrochloric acid)
- highly concentrated alcohols
- chemical denture cleaners
- abrasive cleaning agents

8. Mounting

Bonding extra-orally:

Only use special adhesives, adapted to the material you want to combine with the polycarbonate. Pay attention to the instructions of the adhesive manufacturer.

Cementation:

For temporary cementations temporary cements and for permanent cementations self-adhesive glass ionomer cements are recommended.

9. Disposal

Remains of the polycarbonate blank can be disposed of in the regular household garbage.

We reserve the right to make changes because of the continuous development of our products. Please find the current version of the instruction manual on our website: www.dentaldirekt.de.

This version replaces all previous versions.