Manuel d'instructions

Instruction Manual

DD Bio ZX²

DD Bio ZX²

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DD Bio ZX² 98 (color), DD Bio ZX² 71 (color),
DD Bio ZX² 19, DD Bio ZX² 95, DD Bio ZX² 100

Valid for the variants:

DD Bio ZX² 98 (color), DD Bio ZX² 71 (color)
DD Bio ZX² 19, DD Bio ZX² 95, DD Bio ZX² 100

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Instruction Manual

1. Indications

DD Bio ZX² are dental mill blanks made from zirconium dioxide for the manufacture of dental prostheses. DD Bio ZX² can be used for all dental constructions with up to two bridge elements in the front and side tooth area. Because of the increased translucency, the material is particularly suitable for production of full contour restorations with or without incisal cutback.

2. Contraindications

Para-functional habits (e.g. bruxism), insufficient space, inadequate preparation, known intolerance to the contained components, insufficient oral hygiene.

3. Safety information

Please pay attention to the information in the safety data sheet in its current version.

Avoid the inhalation of dust particles during processing.

Wear protection gloves, safety glasses and a dust mask to prevent irritations of eyes, skin and respiratory system.

4. Handling and storage

The blanks are manufactured and tested to highest quality standards.

Verify the integrity of the packaging and the blank before the first processing. Check if the content of the packaging corresponds to the declaration given on the label. Do not use damaged material.

Storage only in the original containers in a cool and dry environment. Avoid vibrations, contaminations and contact with fluids.

5. Processing / Designing

Handling this medical device should only be allowed to trained staff.

Following construction parameters need to be considered:

Indication	Scheme	Wall thickness for all DD Bio ZX ² variants (excl. ZX ² 19) [mm]	Wall thickness for DD Bio ZX ² Sirona inLab® [mm]	Connector cross section [mm ²]
Single crown	X	0,4	0,5	-
Telescope	X	0,5	0,5	-
Bridge 3-elements				
Anterior	XOX	0,5	0,5	> 7
Posterior		0,5	0,5	> 9
Bridge 4-elements				
Anterior	XOOX	0,6	0,7	> 9
Posterior		0,6	0,7	> 12

Legend: X = abutment tooth O = pontic

Please note: Depending on the construction, the connector cross section might need to be of larger dimensions. An oval connector cross section is desirable; the height of the connector is the decisive factor. Frameworks for ceramic veneering need to be designed in a way to support the veneering ceramic in the cusp region and allow a constant layer thickness. Shoulder or chamfer preparation is recommended.

6. Milling, sintering and further processing

Machining of blanks should always be conducted in the corresponding milling system. All information given by the manufacturer of the machine shall be taken into account. The frameworks need to be examined for visible defects.

Do not process damaged or contaminated frameworks. The white framework can be colored before final sintering with DD Bio ZX² monolith zero (pay attention to separate instruction manual). Pre-colored blanks are not allowed to use with patented coloring technologies.

⚠ Please pay attention to our separate sintering instruction!

Sintering cycle for normal furnace filling without cover:

- ↑ Heating up to 900°C (1652°F) with 8°C/min (46°F/min),
- dwell at 900°C (1652°F) for 30 minutes,
- ↑ heating up to final temp. 1450°C (2642°F) with 3°C/min (37°F/min),
- dwell at final sintering temperature 1450°C (2642°F) for 120 min,
- ↓ cooling to 200°C (392°F) with 10°C/min (50°F/min).

Avoid additional mechanical actions after sintering like blasting or grinding.

If corrections are inevitable use water-cooled tools for the conditioning and make sure that there is no development of heat, which could lead to cracks in the material. Please work only with sharp, diamond-coated tools and very slight contact pressure. Areas that are under tension in clinical use (e.g. connectors) must not be processed. Do not separate in interdental regions. Avoid sharp edges.

7. Mounting

For anchoring we recommend conventional cementation with zinc oxide phosphates cements or glassionomer cements. Also luting composites may be used. Sufficient retentions and a stump with height of at least 3 mm are essential.

We reserve the right to make changes as a result of the continuous development of our products. Please find the current version of the instruction manual on our website: www.dentaldirekt.de This version replaces all previous versions.