

# *DDcubeX<sup>2</sup>*

cubic zirconia system





DD Bio Z classic  
1450°C  
35% translucency



DD Bio ZX² translucent  
1450°C  
40% translucency



DD cubeX² EHT  
1450°C  
49% translucency

# DD cubeX<sup>2</sup>

cubic zirconia system



## The Zirconia (R)evolution

translucent like a lithium disilicate—strength like a zirconia

Did you think there is no more evolution with zirconia? You have been wrong. The new DD cubeX<sup>2</sup> – cubic zirconia system is a material revolution for technicians who love aesthetics and excellence. Zirconia was never like that. A highly biocompatible material with outstanding transparency. The DD cubeX<sup>2</sup> system combines the known positive properties of zirconia with a significant increase in translucency as it is known only from lithium disilicate. This material is perfect for all highly aesthetic treatments in anterior tooth region. Moreover, it is also perfectly appropriate for use in the posterior tooth region. Don't just take our word for it – see for yourself and share our enthusiasm.

